



M I C R O S E N S

Product Overview

MICROSENS MSFET3330 amplifier MSFET-ANALOG

Interface electronics to drive the MSFET3330 pH sensor element



Key Features

- Small packaging
- Sensor status indication
- Switch for reading the integrated temperature diode
- Analogue output

Applications

- Laboratory
- Benchtop setup
- MSFET 3330 Demo Kit

Characteristics

- Temperature range:
-20°C ... 80°C
- Box material:
Polycarbonate

MICROSENS SA

Rue de la Maladière 71c, CP 645, CH-2002 Neuchâtel/Switzerland, Tel: +41-32-724 55 66

EPFL Innovation Park, Batiment D, 1015 Lausanne, Switzerland, Tel.:+41-21-353 91 35

e-mail : info@microsens.ch, www.microsens.ch